

Anti-Human/Mouse/Rat RRN3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05167

3th Edition

Description:RRN3 (RRN3 Homolog, RNA Polymerase I Transcription Factor) is a Protein Coding gene. Among its related pathways are RNA Polymerase I Promoter Escape and mTOR signaling pathway (Pathway Interaction Database). GO annotations related to this gene include binding and RNA polymerase binding.

Antigen:Synthesized peptide derived from human TIF-IA Polyclonal

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: