

Anti-Human/Mouse/Rat MUTYH Polyclonal Antibody
Polyclonal Antibody

Cat.NO.: PA05191

3th Edition

Description:Involved in oxidative DNA damage repair. Initiates repair of A*oxoG to C*G by removing the inappropriately paired adenine base from the DNA backbone. Possesses both adenine and 2-OH-A DNA glycosylase activities.

Antigen:Synthesized peptide derived from the Internal region of human MYH

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: