

Anti-Human/Mouse MTM1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05201

3th Edition

Description: This gene encodes a dual-specificity phosphatase that acts on both phosphotyrosine and phosphoserine. It is required for muscle cell differentiation and mutations in this gene have been identified as being responsible for X-linked myotubular myopathy.

Antigen: Synthesized peptide derived from the Internal region of human Myotubularin.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: