

Anti-Human NOL3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05202

3th Edition

Description: Isoform 1 may be involved in RNA splicing. Isoform 2 may inhibit apoptosis.

Antigen: Synthesized peptide derived from the C-terminal region of human Myp

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Highly expressed in heart and skeletal muscle. Detected at low levels in placenta, liver, kidney and pancreas.