

Anti-Human/Mouse/Rat/Monkey PPP1R12A Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05203

3th Edition

Description:Regulates myosin phosphatase activity.

Antigen:Synthesized peptide derived from human MYPT1 around the non-phosphorylation site of Thr696.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expressed in striated muscles, specifically in type 2a fibers (at protein level).