

Anti-Human/Mouse/Rat BCL6 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05207

3th Edition

Description: Transcriptional repressor which is required for germinal center formation and antibody affinity maturation. Probably plays an important role in lymphomagenesis.

Antigen: Synthesized peptide derived from B-cell lymphoma 6 protein

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in germinal center T and B cells and in primary immature dendritic cells.