

Anti-Human N4BP1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05212

3th Edition

Description:N4BP1 (NEDD4 Binding Protein 1) is a Protein Coding gene. GO annotations related to this gene include RNA binding. An important paralog of this gene is KHNYN.

Antigen:Synthesized peptide derived from the Internal region of human N4BP1

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Detected in heart, lung, brain, liver, skeletal muscle, pancreas, kidney, spleen, testis and ovary.