

**Anti-Human/Mouse/Rat Cleaved-F2R (S42) Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05216**

---

3th Edition

**Description:**Coagulation factor II receptor is a 7-transmembrane receptor involved in the regulation of thrombotic response. Proteolytic cleavage leads to the activation of the receptor. F2R is a G-protein coupled receptor family member.

**Antigen:**Synthesized peptide derived from human PAR1 (Cleaved-Ser42) Polyclonal

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Platelets and vascular endothelial cells.