

Anti-Human/Mouse/Rat TMEM72 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05219

3th Edition

Description: This gene encodes a transmembrane protein which may be expressed specifically in the kidney. TMEM72 (Transmembrane Protein 72) is a Protein Coding gene.

Antigen: Synthesized peptide derived from human Tmem72 Polyclonal

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: