

**Anti-Human/Mouse NDUFA4L2 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05236**

---

3th Edition

**Description:**NDUFA4L2 (NDUFA4, Mitochondrial Complex Associated Like 2) is a Protein Coding gene. Among its related pathways are Respiratory electron transport, ATP synthesis by chemiosmotic coupling, and heat production by uncoupling proteins. and GABAergic synapse. An important paralog of this gene is NDUFA4.

**Antigen:**Synthesized peptide derived from the C-terminal region of human NDUFA4L2

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**