

**Anti-Human/Mouse/Rat NEGR1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05245**

---

3th Edition

**Description:**NEGR1 (Neuronal Growth Regulator 1) is a Protein Coding gene. Diseases associated with NEGR1 include Podoconiosis. Among its related pathways are Post-translational modification- synthesis of GPI-anchored proteins and Metabolism of proteins. An important paralog of this gene is LSAMP.

**Antigen:**Synthesized peptide derived from the Internal region of human NEGR1

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**