

**Anti-Human/Mouse SYNE3 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05250**

---

3th Edition

**Description:**SYNE3 (Spectrin Repeat Containing Nuclear Envelope Family Member 3) is a Protein Coding gene. GO annotations related to this gene include actin binding and actin filament binding.

**Antigen:**Synthesized peptide derived from the Internal region of human Nesprin-3

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Expressed in aortic endothelial cells (at protein level).