

Anti-Human/Mouse/Rat NRSN1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05260

3th Edition

Description:NRSN1 (Neurensin 1) is a Protein Coding gene. Diseases associated with NRSN1 include Louse-Borne Relapsing Fever. An important paralog of this gene is NRSN2.

Antigen:Synthesized peptide derived from the Internal region of human Neurensin-1

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expressed in brain. Not detectable in other tissues tested.