

Anti-Human/Rat NXPH4 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05262

3th Edition

Description: NXPH4 (Neurexophilin 4) is a Protein Coding gene. Diseases associated with NXPH4 include Joubert Syndrome 1. An important paralog of this gene is NXPH2.

Antigen: Synthesized peptide derived from the C-terminal region of human Neurexophilin-4

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in brain, spleen, and testis.