

Anti-Human ANO7 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05279

3th Edition

Description: Does not exhibit calcium-activated chloride channel (CaCC) activity. May play a role in cell-cell interactions.

Antigen: Synthesized peptide derived from the C-terminal region of human NGEP

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Specifically expressed in epithelial cells of the prostate (at protein level).