



Anti-Human/Mouse/Rat NGF Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05280

3th Edition

Description: Nerve growth factor is important for the development and maintenance of the sympathetic and sensory nervous systems. It stimulates division and differentiation of sympathetic and embryonic sensory neurons.

Antigen: Synthesized peptide derived from the Internal region of human NGF.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: