

Anti-Human/Mouse NUP160 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05328

3th Edition

Description: NUP160 is 1 of up to 60 proteins that make up the 120-MD nuclear pore complex, which mediates nucleoplasmic transport. NUP160 (Nucleoporin 160) is a Protein Coding gene. Among its related pathways are HIV Life Cycle and Transport of the SLBP independent Mature mRNA. GO annotations related to this gene include nucleocytoplasmic transporter activity.

Antigen: Synthesized peptide derived from the Internal region of human Nup160

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: