

Anti-Human/Mouse/Rat Oct-1/2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05339

3th Edition

Description:The OCT1 transcription factor was among the first identified members of the POU transcription factor family (summarized by Sturm et al., 1993 [PubMed 8314572]). Members of this family contain the POU domain, a 160-amino acid region necessary for DNA binding to the octameric sequence ATGCAAAT

Antigen:Synthesized peptide derived from human 41183 around the non-phosphorylation site of S385.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: