

**Anti-Human/Mouse Oct-1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05340**

---

3th Edition

**Description:**The OCT1 transcription factor was among the first identified members of the POU transcription factor family. Members of this family contain the POU domain, a 160-amino acid region necessary for DNA binding to the octameric sequence ATGCAAAT.

**Antigen:**Synthesized peptide derived from the N-terminal region of human Oct-1

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Ubiquitous. Isoform 2 is lymphocyte-specific.