

Anti-Human Olfactory receptor 10V1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05361

3th Edition

Description:Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. OR10V1 (Olfactory Receptor Family 10 Subfamily V Member 1) is a Protein Coding gene. Among its related pathways are Olfactory Signaling Pathway and Signaling by GPCR. GO annotations related to this gene include G-protein coupled receptor activity and odorant binding. An important paralog of this gene is OR10Q1.

Antigen:Synthesized peptide derived from the Internal region of human Olfactory receptor 10V1

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: