

Anti-Human Olfactory receptor 2AG1/2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05376

3th Edition

Description:OR2AG1 (Olfactory Receptor Family 2 Subfamily AG Member 1 (Gene/Pseudogene)) is a Protein Coding gene. Among its related pathways are Peptide ligand-binding receptors and Olfactory Signaling Pathway. GO annotations related to this gene include G-protein coupled receptor activity and olfactory receptor activity. An important paralog of this gene is OR2AG2. OR2AG2 (Olfactory Receptor Family 2 Subfamily AG Member 2) is a Protein Coding gene. Among its related pathways are Olfactory Signaling Pathway and Signaling by GPCR. GO annotations related to this gene include G-protein coupled receptor activity and olfactory receptor activity. An important paralog of this gene is OR2AG1.

Antigen:Synthesized peptide derived from the Internal region of human Olfactory receptor 2AG1/2

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: