

Anti-Human Olfactory receptor 2T1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05384

3th Edition

Description:OR2T1 is an olfactory receptor. Members of this family interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) which arise from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for humans are independent to those assigned to other organisms.

Antigen:Synthesized peptide derived from the C-terminal region of human Olfactory receptor 2T1

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: