

Anti-Human/Mouse/Rat p140 Cap Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05468

3th Edition

Description: SRCIN1 (SRC Kinase Signaling Inhibitor 1) is a Protein Coding gene. Diseases associated with SRCIN1 include Sclerosing Keratitis and Punctate Epithelial Keratoconjunctivitis. Among its related pathways are Adhesion and Neuroscience. GO annotations related to this gene include protein kinase binding and protein domain specific binding. An important paralog of this gene is KIAA1217.

Antigen: Synthesized peptide derived from the Internal region of human p140Cap

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in some primary breast carcinomas where its presence is significantly associated with increased tumor size. Not detected in normal breast tissue.