

Anti-Human/Mouse/Rat p15 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05469

3th Edition

Description: Interacts strongly with CDK4 and CDK6. Potent inhibitor. Potential effector of TGF-beta induced cell cycle arrest.

Antigen: Synthesized peptide derived from the C-terminal region of human p15

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Isoform 2 is expressed in normal (keratinocytes, fibroblasts) and tumor cell lines.