

**Anti-Human/Mouse/Monkey p18 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05471**

---

3th Edition

**Description:** Interacts strongly with CDK6, weakly with CDK4. Inhibits cell growth and proliferation with a correlated dependence on endogenous retinoblastoma protein RB.

**Antigen:** Synthesized peptide derived from the C-terminal region of human p18

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Highest levels found in skeletal muscle. Also found in pancreas and heart.