

Anti-Human/Mouse/Rat p35 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05482

3th Edition

Description:p35 is a neuron specific activator of CDK5. The complex p35/CDK5 is required for neurite outgrowth and cortical lamination. Activator of TPKII.

Antigen:Synthesized peptide derived from the N-terminal region of human p35

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Brain and neuron specific.