

Anti-Human/Mouse/Rat PAH Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05507

3th Edition

Description:PAH encodes the enzyme phenylalanine hydroxylase that is the rate-limiting step in phenylalanine catabolism. Deficiency of this enzyme activity results in the autosomal recessive disorder phenylketonuria.

Antigen:Synthesized peptide derived from the C-terminal region of human PAH.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: