

Anti-Human/Mouse/Rat SERPINE1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05508

3th Edition

Description: Serine protease inhibitor. This inhibitor acts as 'bait' for tissue plasminogen activator, urokinase, protein C and matrilysin-3/TMPRSS7. Its rapid interaction with PLAT may function as a major control point in the regulation of fibrinolysis.

Antigen: Synthesized peptide derived from the Internal region of human PAI-1.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Found in plasma and platelets and in endothelial, hepatoma and fibrosarcoma cells.