

Anti-Human/Mouse PAK5/6 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05511

3th Edition

Description:The activated kinase acts on a variety of targets.

Antigen:Synthesized peptide derived from human PAK5/6 around the non-phosphorylation site of S602/S560.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:PAK5 is Predominantly expressed in brain.PAK6 isSelectively expressed in brain and testis, with lower levels in multiple tissues including prostate and breast.