

**Anti-Human/Mouse/Rat/Monkey Paxillin Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05523**

---

3th Edition

**Description:** Cytoskeletal protein involved in actin-membrane attachment at sites of cell adhesion to the extracellular matrix (focal adhesion).

**Antigen:** Synthesized peptide derived from human Paxillin around the non-phosphorylation site of Tyr31.

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**