

**Anti-Human/Mouse PDIK1L Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05538**

---

3th Edition

**Description:**PDIK1L (PDLIM1 Interacting Kinase 1 Like) is a Protein Coding gene. GO annotations related to this gene include transferase activity, transferring phosphorus-containing groups and protein tyrosine kinase activity. An important paralog of this gene is STK35.

**Antigen:**Synthesized peptide derived from the Internal region of human PDIK1L

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Expressed in liver, kidney, pancreas, spleen, thymus and prostate.