

Anti-Human/Mouse/Rat PDRG1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05543

3th Edition

Description:PDRG1 (P53 And DNA Damage Regulated 1) is a Protein Coding gene. GO annotations related to this gene include unfolded protein binding.

Antigen:Synthesized peptide derived from the Internal region of human PDRG1

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Predominantly expressed in normal testis and exhibits reduced but detectable expression in other organs.