

Anti-Human/Mouse/Rat PDXDC1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05544

3th Edition

Description:PDXDC1 (Pyridoxal Dependent Decarboxylase Domain Containing 1) is a Protein Coding gene. GO annotations related to this gene include pyridoxal phosphate binding and carboxy-lyase activity.

Antigen:Synthesized peptide derived from the C-terminal region of human PDXDC1

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: