

Anti-Human/Mouse N4BP2L2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05557

3th Edition

Description:N4BP2L2 (NEDD4 Binding Protein 2 Like 2) is a Protein Coding gene. GO annotations related to this gene include enzyme binding and RNA polymerase II transcription corepressor activity. An important paralog of this gene is N4BP2.

Antigen:Synthesized peptide derived from the Internal region of human PFAAP5

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: