

Anti-Human PGLS Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05567

3th Edition

Description:PGLS (6-Phosphogluconolactonase) is a Protein Coding gene. Diseases associated with PGLS include Lymphatic System Disease. Among its related pathways are Pentose phosphate pathway and Pentose Phosphate Pathway (Erythrocyte). GO annotations related to this gene include monosaccharide binding and 6-phosphogluconolactonase activity.

Antigen:Synthesized peptide derived from the C-terminal region of human PGLS

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: