

Anti-Human/Mouse/Rat PHLDA3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05573

3th Edition

Description:PHLDA3 (Pleckstrin Homology Like Domain Family A Member 3) is a Protein Coding gene. Among its related pathways are Apoptosis and Autophagy and PI3K / Akt Signaling. GO annotations related to this gene include phosphatidylinositol-4,5-bisphosphate binding and phosphatidylinositol-3-phosphate binding. An important paralog of this gene is PHLDA1.

Antigen:Synthesized peptide derived from the Internal region of human PHLDA3

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Widely expressed with lowest expression in liver and spleen.