

Anti-Human OLFML2A Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05574

3th Edition

Description:OLFML2A (Olfactomedin Like 2A) is a Protein Coding gene. GO annotations related to this gene include protein homodimerization activity and extracellular matrix binding. An important paralog of this gene is OLFML2B.

Antigen:Synthesized peptide derived from the Internal region of human Photomedin-1

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:In the kidney expressed only by podocytes, wherein they localize to major processes.