

Anti-Human/Mouse OLFML2B Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05575

3th Edition

Description: This gene encodes an olfactomedin domain-containing protein. Most olfactomedin domain-containing proteins are secreted glycoproteins. OLFML2B (Olfactomedin Like 2B) is a Protein Coding gene. GO annotations related to this gene include extracellular matrix binding. An important paralog of this gene is OLFML2A.

Antigen: Synthesized peptide derived from the Internal region of human Photomedin-2

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: