

Anti-Human PIPOX Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05588

3th Edition

Description:PIPOX (Pipecolic Acid And Sarcosine Oxidase) is a Protein Coding gene. Diseases associated with PIPOX include Peroxisomal Biogenesis Disorders. Among its related pathways are Lysine degradation and Peroxisome. GO annotations related to this gene include receptor binding and sarcosine oxidase activity.

Antigen:Synthesized peptide derived from the Internal region of human PIPOX

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: