

Anti-Human PJA2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05591

3th Edition

Description:PJA2 (Prajna Ring Finger Ubiquitin Ligase 2) is a Protein Coding gene. Among its related pathways are Innate Immune System and Class I MHC mediated antigen processing and presentation. GO annotations related to this gene include ligase activity and protein kinase A catalytic subunit binding. An important paralog of this gene is PJA1.

Antigen:Synthesized peptide derived from the Internal region of human PJA2

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: