

**Anti-Human PMEPA1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05603**

---

3th Edition

**Description:**Involved in down-regulation of the androgen receptor (AR). Enhances ubiquitination and proteasome-mediated degradation of AR, probably by recruiting NEDD4.

**Antigen:**Synthesized peptide derived from the Internal region of human PMEPA1

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Highest expression in prostate. Also expressed in ovary.