

Anti-Human/Monkey PXMP4 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05605

3th Edition

Description:PXMP4 (Peroxisomal Membrane Protein 4) is a Protein Coding gene. Diseases associated with PXMP4 include Prostate Adenocarcinoma. Among its related pathways are Peroxisome.

Antigen:Synthesized peptide derived from the N-terminal region of human Pmp24

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expressed in normal prostate epithelial cells, and androgen-sensitive prostate adenocarcinoma cells. Not expressed in androgen-insensitive prostate adenocarcinoma cells.