

**Anti-Human/Mouse PNPT1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05608**

---

3th Edition

**Description:**PNPT1 is a subunit of the exosome complex, which is involved in 3-prime-to-5-prime exoribonuclease activity for RNA processing and degradation. It contains one KH domain and one S1 motif domain.

**Antigen:**Synthesized peptide derived from the C-terminal region of human PNPase

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**