

**Anti-Human/Mouse/Rat POLR3H Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05612**

---

3th Edition

**Description:** POLR3H (RNA Polymerase III Subunit H) is a Protein Coding gene. Among its related pathways are Innate Immune System and RNA Polymerase III Transcription Initiation. GO annotations related to this gene include DNA-directed 5-3 RNA polymerase activity and RNA polymerase III activity.

**Antigen:** Synthesized peptide derived from the C-terminal region of human POLR3H

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**