

Anti-Human/Monkey PPP4R1L Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05620

3th Edition

Description:PPP4R1L (Protein Phosphatase 4 Regulatory Subunit 1 Like (Pseudogene)) is a Pseudogene.

Antigen:Synthesized peptide derived from the Internal region of human PPP4R1L

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: