

Anti-Human/Mouse/Rat PRPF19 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05631

3th Edition

Description: PSO4 is the human homolog of yeast Pso4, a gene essential for cell survival and DNA repair (Beck et al., 2008 [PubMed 18263876]).

Antigen: Synthesized peptide derived from the Internal region of human PRP19

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Ubiquitous. Weakly expressed in senescent cells of different tissue origins. Highly expressed in tumor cell lines.