

Anti-Human/Mouse/Rat PSMD12 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05643

3th Edition

Description: Acts as a regulatory subunit of the 26S proteasome which is involved in the ATP-dependent degradation of ubiquitinated proteins.

Antigen: Synthesized peptide derived from the Internal region of human PSMD12

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: