

Anti-Human RAB40B Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05671

3th Edition

Description:RAB40B is similar to a yeast protein which suggests a role of the gene product in regulating secretory vesicles. It may be a substrate recognition component of a SCF like ECS (Elongin Cullin SOCS box protein) E3 ubiquitin ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins.

Antigen:Synthesized peptide derived from the C-terminal region of human Rab 40B

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: