

Anti-Human/Mouse/Rat RAB6A Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05672

3th Edition

Description:Protein transport. Regulator of membrane traffic from the Golgi apparatus towards the endoplasmic reticulum (ER). Has a low GTPase activity.

Antigen:Synthesized peptide derived from the Internal region of human Rab 6A

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Ubiquitous.