

**Anti-Human/Mouse RAB11FIP4 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05677**

---

3th Edition

**Description:** Proteins of the large Rab GTPase family (see RAB1A; MIM 179508) have regulatory roles in the formation, targeting, and fusion of intracellular transport vesicles. RAB11FIP4 is one of many proteins that interact with and regulate Rab GTPases (Hales et al., 2001 [PubMed 11495908]). The RAB11FIP4 effector protein acts in endosomal trafficking.

**Antigen:** Synthesized peptide derived from the Internal region of human Rab11-FIP4

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Present at high level in testis (at protein level). Weakly expressed in other tissues.